

394Δ 8710Δ

Brewer Spectrophotometer laboratory wavelength characterization with a tunable-laser implications to the Brewer ozone retrieval

Alberto Redondas^{1,2}, Saulius Nevas³, Alberto Berjón^{4,2}, Meelis-Mait Sildoja³, and Sergio F. León-Luis^{1,2}

¹Agencia Estatal de Meteorología, Izaña Atmospheric Research Center, Spain

²Regional Brewer Calibration Center for Europe, Izaña Atmospheric Research Center, Tenerife, Spain

³Physikalisch-Technische Bundesanstalt (PTB), Braunschweig and Berlin, Germany

⁴University of La Laguna, Department of Industrial Engineering, S.C. de Tenerife, Spain

Correspondence to: Alberto Redondas (aredondasm@aemet.es)

Abstract.

In this work we present the wavelength calibration of the traveling reference Brewer of the RBCC-E <http://rbcce.aemet.es> at PTB in Braunschweig, Germany. The wavelength calibration is needed for calculation of the ozone absorption coefficient used by the Brewer ozone algorithm. In order to validate the standard procedure for determining Brewers wavelength calibration an experiment using an tunable laser light source have been performed within the ATMOZ project. In this experiment we compare the standard procedure of wavelength calibration performed to Brewer instrument uses to the calibration results obtained using the PTB laser facility. This allows to check the methodology measuring directly the ozone absorption coefficient without the assumptions made by the operational methodology. The results of the laser experiment reproduces the the operational methodology used and shows that there is a underestimation of 0.8% due the use of the parametrized slit function.

half maximum) and the relation between the positions of the grating and the corresponding wavelengths (dispersion relation) for the operational wavelengths used for the ozone determination. The Brewer spectrophotometer have two operating modes. The ozone one, used for ozone measurement, is performed with the diffraction grating at a fixed position while the six operational wavelengths are measured by rotating slit mask. The scanning mode, used for the spectral UV measurements, is performed while the slits are fixed and the spectral scan is carried out by turning the diffraction grating. To obtain the ozone absorption coefficient the instrumental slit function is convoluted with the Bass & Paur ozone cross section. The use of the of the laser tunable source allow us to:

- Calculate the ozone absorption coefficient directly. The ozone absorption coefficient determination uses the scan the spectral lines so that dispersion relation is used to convert grating position in micrometer steps to wavelengths assuming a quadratic relation. Scanning with the laser around the ozone operational wavelengths we can determine the instrumental slit functions directly and weight them with the ozone cross sections without need for the assumptions of the slit functions and the dispersion relations used in the normal operational procedure.
- Calculate the dispersion relation based on regularly spaced reference spectral lines provided by the tunable laser instead of the irregularly distributed emission lines of the Hg, Cd and Zn spectral lamps.

During the experiment we perform three measurements:

Copyright statement. TEXT

1 Background

The wavelength calibration is needed for calculation of the ozone absorption coefficient used by the Brewer ozone algorithm. The usual wavelength calibration procedure is performed by analyzing recorded emission lines of the spectral discharge lamps usually Mercury (Hg), Cadmium (Cd), and Zinc (Zn) detailed on Table 1, determining their central wavelengths and the corresponding FWHM (full width

1. The standard method of the dispersion measurements using spectral lamps described on section 2.
2. Direct dispersion measurements (Laser scanning). While Brewer is measuring in ozone mode and in aerosol mode, the laser will scan around (± 2 nm) the six Brewer slits with a step of 0.05nm for different grating positions
3. Dispersion measurements using tunable laser (Brewer scanning) : While Laser is emitting in a fixed wavelength, the Brewer will scan around this wavelength (± 2 nm) moving the grating and using the 6 slits. This process will be repeated 16 times, using wavelengths ranging from 290 nm to 365 nm with an increment of 5nm. This allow us to estimate the error due to the lack of spectral lines on the end wavelength range of the brewer spectrometer and the fact that that the emission lines are not equally spaced.

2 The Brewer spectrometer calibration

The Brewer instrument measures the intensity of direct sunlight at six wavelengths (λ) in the UV (303.2, 306.3, 310.1, 313.5, 316.8, and 320.1 nm) each covering a bandwidth of 0.5 nm (resolution power $\lambda/\delta\lambda$ of around 600). The spectral measurement is achieved by a holographic grating in combination with a slit mask which selects the channel to be analyzed by a photomultiplier. The longest four wavelengths are used for the ozone calculation. Based on the Lambert-Beer law, the total ozone column in the Brewer algorithm can be expressed as:

$$X = \frac{F - ETC}{\alpha\mu} \quad (1)$$

where F are the measured double ratios corrected for Rayleigh effects, α is the ozone absorption coefficient, μ is the ozone air mass factor, and ETC is the extra-terrestrial constant. The F , α and ETC parameters are weighted functions at the operational wavelengths with weighting coefficients w :

$$F = \sum_i^4 w_i F_i - \frac{p}{p_0} \beta_i \mu \quad (2)$$

$$\alpha = \sum_i^4 w_i \alpha_i \quad (3)$$

$$ETC = \sum_i^4 w_i F_{0i} \quad (4)$$

where, β_i are the Rayleigh coefficients, p is the climatological pressure at the measurement site, p_0 is the pressure at sea level, and F_0 are the individual extra-terrestrial constants at each wavelength. The weights $w = [1, -0.5, -2.2, 1.7]$ have been chosen so as to minimize the influence of SO_2 and verify:

$$\sum_i^4 w_i = 0 \quad (5)$$

$$\sum_i^4 w_i \lambda_i = 0 \quad (6)$$

This widely eliminates absorption features which depend, in local approximation, linearly on the wavelength, like for example the contribution from aerosols.

We can divide the calibration in three steps:

1. Instrumental, wavelength, and ETC transfer: the Instrumental calibration includes all the parameters that affect the measured counts (F), in particular Dead Time correction, Temperature coefficients and Filter attenuation.
2. Wavelength calibration to determine the ozone absorption coefficient: the so-called “dispersion test” are used to obtain the particular wavelength for the instrument and the slit, or instrumental function, of each spectrophotometer. Note that the precise wavelengths of every Brewer spectrophotometer are slightly different from instrument to instrument.
3. Finally, the ETC transfer is performed by comparison with the reference or, in the case of the reference instruments, by the Langley method.

The Brewer wavelength calibration follows the operative procedure (Gröbner et al., 1998; Kerr, 2002) used by the Regional Brewer Calibration Center-Europe (RBCC-E) at the calibration campaigns. In summary the individual wavelengths in the Brewer measurement set are located at the focal plane of the spectrometer through the use of a stainless steel mask of seven slits, the particular wavelength is determined by analyzing the measurements of a series of discharge lamps on a process so-called dispersion-test which determines the central wavelength and Full Width Half Maximum (FWHM) of every slits. Then the wavelength setting is optimized to minimize the effect of wavelength shift during the operation of the instrument by the sun-scan test (Lamb, K and Asbridge, A.I. , 1996). And finally the ozone absorption coefficient is determined.

The ozone absorption coefficient is defined as:

$$\tilde{\alpha}(X, \mu) = \sum w_i \frac{\int \alpha(\lambda) * S(\lambda, \lambda') * F(\lambda, \lambda', X, \mu) d\lambda'}{\int S(\lambda, \lambda') * F(\lambda, \lambda', X, \mu) d\lambda'} \quad (7)$$

Where S is the instrument slit function for the corresponding wavelength F is the sun spectra mostly which depends on the UV. range of the ozone concentration and airmass, and σ the ozone cross section at temperature, -46.3°C for Dobson

5 Network and -45°C for Brewer instruments.

The Brewer operative method, perform the following assumptions.

1. Use ideal slits, the slits function are parametrized as a trapezoid; and isosceles triangle truncated at 0.82 height.
2. Stray light is not considered the transmission is outside the triangle.
3. The FWHM of the triangle is dependent of the slit as is derived from the dispersion test.
- 15 4. The ozone cross section is Bass & Paur absorption coefficient.
5. Solar Spectra is not considered , ($F=1$)

With this assumption the ozone effective absorption is essentially the same as the approximation method of Bernhard

20 et al. (2005) used on Dobson Spectrometer. (see Eq. 8).

$$\alpha_i = \frac{\int \sigma(\lambda) S_i(\lambda) d\lambda}{\int S_i(\lambda) d\lambda} \quad (8)$$

2.1 Dispersion Test

The Brewer spectrophotometer is constructed with either one or two modified EbertFastie type monochromators single or

25 double Brewer. The first monochromator disperses the incoming radiation onto six exit slits. In the case of the double Brewer , the six exit slits (intermediate slits) of the first monochromator are the entrance slits to a second monochromator that is used in the recombining mode. The wavelength

30 is selected by choosing one of the six exit slit (ozone mode) or rotating the grating (uv mode). The rotation of the grating is managed by a micrometer. The drive mechanism, consisting of a motor-driven micrometer linked to an arm that rotates the grating. The smallest wavelength increment

35 one motor step varies steadily from approximately 8.0 pm to 7.0 pm (0.0080 nm).

The dispersion relation (which gives the relation between the micrometer step and wavelength) is determined scanning the emission lines. The brewer scan every 10 steps (0.6A),

40 the central step and the FWHM are calculated assuming a isosceles triangle. The both sides of the peak are fitted to straight line ,taking only the observations above 20% and below 80% of the normalized peak, the central point is calculated by the intersection point and the FWHM is the with

45 of the triangle. (Figure x). Then the dispersion relation is calculated as quadratic polynomial or using the cubic approximation (Gröbner et al., 1998). This relation is used to trans-

form the previously determined central and FWHM in micrometer steps to wavelength scale.

The Stability of the wavelength calibration during the brewer operation is checked by the measure of the internal HG line, in most of the brewer a 302 nm double line is used, (302.150 nm and 302.347 nm) due to its proximity to the operational wavelengths however on brewer #185 and in increased number of brewer the test is performed using the more powerful 296.7 nm line. The test taking 12 measurements, 10 steps of the micrometer motor apart, of the intensity from the mercury lamp on slit 0. The registered peak is compared to a reference stored one. The comparison the scans is done by shifting the two scans against each other and calculating the correlation coefficient between the two after each shift. The interpolated step number resulting in the maximum of the correlation is the reference micrometer position. (Savastiouk, 2005) . If the required adjustment of the micrometer position is more than one and a half steps, the test is repeated. Due to the hg test the accuracy of the wavelength setting cannot be better than about micrometer step (0.86) which, determined by is limited by one and half step correction, this translates in ozone absorption coefficient change approximately of 0.10 atm cm^{-1} which gives approximately 0.3% in ozone concentration.

2.2 Sun Scan

Savastiouk (2005)

3 OPO experiment

3.1 Laser Instrumental Setup

The facility PLACOS (Pulsed Laser for Advanced Characterisation of Spectroradiometers) at Physikalisch-Technische Bundesanstalt (PTB). This system is employed for the spectral characterisation of array spectroradiometers throughout the spectral range from 220 nm to 2200 nm. It uses an optical parametric oscillator (OPO) and a second-harmonic generator (SHG). The pump radiation for the OPO crystals is delivered by the 3rd harmonic of a Nd:YAG laser (355 nm). The resulting signal of the OPO covers the spectral range from approximately 410 nm up to 2200 nm. The SHG unit enables an extended range down to 220 nm. (Figure 1)

In contrast with the standard procedure where the Brewer scans the spectral lamps lines, on this experiment the brewer measure in ozone mode, the grating is fixed on the ozone position and measure the laser system with the seven slits. The opo scan every 0.01 nm This experiment measure the seven slits (the slit #1 is blocked) were the OPO system . The experiment is complemented using the laser as source of spectral lines covering the range from 290 nm to 360 nm in regular steps of 5 nm, the brewer resolution of the scan is doubled on the experiment and is done every 5 nm.

3.1.1 PMT No linearity

The brewer photomultiplier (PMT) is not linear with pulsed sources, the PMT manual advised we must change PMT electronics configuration for pulsed sources. As the main objective is to validate the operational wavelength calibration of the brewer we decide to maintain the instrument configuration as operated in field and correct the no-linearity using the simultaneous measures of the laser intensity.

The ratios between the measured counts of the brewer to the recorded laser measurements are shown on figure 2. The no linearity is evident with regions with hysteresis, near 10^4 Brewer counts/seconds, regions where the correction is not reliable around 10^3 and with counts lower than 100. As we can control the intensity of the laser, we are able to work on the flats regions and apply the determined the nonlinear correction. This correction will not affect to the calculation of the central wavelength but affect the determination of FWHM (Figure 3) if not correction is applied.

One of the limitations of this correction is to work with low signals, the correction is not reliable for counts lower than 100. Also we observe (Figure 5) as the recorder dark counts (measurements performed with the blocked slit #1), were highest immediately after exposing the PMT to the laser light. The dark signal was then gradually decreasing with time. The dark values of the PMT may show a fast decay after the excitation, which may cause the values for slit #1 (measured immediately after slit #0) be higher than for some other slits, measured later.

3.1.2 Slit parametrization

The brewer algorithms assume a trapezoidal slit, cut a 0.87 of the height (Figure 4), with the center and fwhm calculated for every slit by the dispersion. This experiment allow us to estimate the difference on the ozone calculation if we use the measured slit instead of the parametrised one. We calculate the ozone absorption coefficient for the four ozone cross section evaluated in the "ACSO" comitee ("Absorption Cross Sections of Ozone") (Orphal et al. (2016).

Two of this are versions of Bass and Paur (1985) cross-sections denoted as Brewer operational (Brw), IGACO quadratic coefficient (B&P), the xs of Daumont Brion Malicet (DBM) (Daumont et al. (1992), Brion et al. (1993), and Malicet et al. (1995)), and the newly recommended data set for ozone ground based calculation Serdyuchenko, Groshelev, Weber (SGW) (Serdyuchenko et al., 2012; Groshelev et al., 2012; Weber et al., 2016).

Using the measured slit, the effective ozone cross section calculated is 0.9% higher respect the parametrized brewer slit on the standard procedure (Table 2) independelty of the cross section used.

4 Comparison

The experiment allow us to validate the Brewer standard methodology to perform the wavelength calibration. Using as a reference the direct calculation we compare the operative method using the spectral lamps and using the laser in regular intervals in both methods using the quadratic fit and the cubic fit.

The figure 7 shows the difference on the central wavelength calculation the quadratic methods shows discrepancies bigger than 0.1 Å for wavelength adobe 320 nm and much bigger near 350 nm. This is also indicated by the residuals of the fit quadratic systematic and much bigger than the cubic method (not shown) that shows the quadratic method is only valid on the ozone range (310-320 nm). The comparison of the FWHM calculation (Figure 7) show a different pattern with a 0.1Å difference on the ozone range and with less difference between quadratic and cubic methods.

The results of the difference to the direct calculation of the ozone absorption coefficient are summarized on the table x, for the six 3

5 Conclusions

1. Using the measured slit instead of the parametrized slit increase the ozone absorption coefficient and consequently the calculated ozone in 0.8%.
2. The quadratic fit used on the standard algorithm is not suitable outside the ozone range 310-320 nm and shows systematic residuals this is particularly important on the higher range
3. The comparison of three experiments shows a maximum difference of 0.3% if the cubic fit is used, the error in the ozone using direct measurement in ozone mode and operative discharge lamp method is only of 0.1%. This confirms the standard procedure using on the RBCC-E calibrations.

Code availability. TEXT

Data availability. TEXT

Code and data availability. TEXT

Acknowledgements. This work has been supported by the European Metrology Research Programme (EMRP) within the joint research project ENV59 Traceability for atmospheric total column ozone (ATMOZ). The EMRP is jointly funded by the EMRP participating countries within EURAMET and the European Union.




Figure 1. TEXT




Figure 2. Log-log plot of the normalized ratio of Brewer registered counts to the Laser counts against the brewer counts, the black points shows the experimentally points and the red curve the fit used to correct the brewer measurements.

References

- Bernhard, G., Evans, R. D., Labow, G. J., and Oltmans, S. J.: Bias in Dobson total ozone measurements at high latitudes due to approximations in calculations of ozone absorption coefficients and air mass, *Journal of Geophysical Research: Atmospheres*, 110, n/a–n/a, <https://doi.org/10.1029/2004JD005559>, <http://dx.doi.org/10.1029/2004JD005559>, d10305, 2005.
- Brion, J., Chakir, A., Daumont, D., Malicet, J., and Parisse, C.: High-resolution laboratory absorption cross section of O₃. Temperature effect, *Chem. Phys. Lett.*, 213, 610–612, [https://doi.org/10.1016/0009-2614\(93\)89169-I](https://doi.org/10.1016/0009-2614(93)89169-I), <http://www.sciencedirect.com/science/article/pii/000926149389169I>, 1993.




Figure 3. Measurement of the 310 slit (Slit #4) uncorrected counts per second (black) and no linearly corrected counts, the central wavelength determination is not affected but FWHM is bigger if the correction is not applied, the hysteresis is evident on the asymmetry of the low intensity of the plot, as the method use only the observations between 0.2 and 0.8 of the normalized intensity do not have effect on the calculations.




Figure 4. Plot of the parametrized slits and the measured slits (dots) and the different ozone cross sections (left axis) used on effective ozone absorption coefficient brewer ozone calculation.

- Daumont, Brion, J., Charbonnier, J., and Malicet, J.: Ozone UV spectroscopy I: Absorption cross-sections at room temperature, *Journal of Atmospheric Chemistry*, 15, 145–155, <https://doi.org/10.1007/BF00053756>, 1992.
- Gorshchev, V., Serdyuchenko, A., Weber, M., and Burrows, J. P.: High spectral resolution ozone absorption cross-sections: Part I. Measurements, data analysis and comparison around 293K., *J. Quant. Spectrosc. Radiat. Transfer*, http://www.iup.uni-bremen.de/UVSAT_material/manuscripts/serdyuchenko_o3x_part1_v1.1.pdf, 2012.
- Gröbner, J., Wardle, D. I., McElroy, C. T., and Kerr, J. B.: Investigation of the wavelength accuracy of




Figure 5. Measurement of the Brewer in ozone mode while the Laser is scanning every 0.1 nm, in blue the measurements of the slit 1 correspond with the dark count.




Figure 6. Measurement of the Brewer in ozone mode while the Laser is scanning every 0.1 nm, in blue the measurements of the slit 1 correspond with the dark count.

Brewer spectrophotometers, Appl. Opt., 37, 8352–8360, <https://doi.org/10.1364/AO.37.008352>, 1998.

Kerr, J.: New methodology for deriving total ozone and other atmospheric variables from Brewer spectrophotometer direct sun spectra, Journal of Geophysical Research: Atmospheres, 107, 2002.

Lamb, K and Asbridge, A.I. : Brewer Technical Note No.96.01 Sun scans and relationship to the wavelength calibration step number, Tech. rep., International Ozone Services, 1996.

Malicet, Daumont, D., Charbonnier, J., Parris, C., Chakir, A., and Brion, J.: Ozone UV spectroscopy. II. Absorption cross-sections and temperature dependence, J. Atmos. Chem, 21, 263–273, <https://doi.org/10.1007/BF00696758>, 1995.




Figure 7. Central Wavelength determination differences between the direct determination and the scanning methods: with the laser in equally spaced lines every 5nm and the discharge lamps using in both cases quadratic and cubic fitting




Figure 8. FWHM determination differences between the direct determination and the scanning methods: with the laser in equally spaced lines every 5nm and the discharge lamps using in both cases quadratic and cubic fitting

Orphal, J., Staehelin, J., Tamminen, J., Braathen, G., De Backer, M.-R., Bais, A., Balis, D., Barbe, A., Bhartia, P. K., Birk, M., et al.: Absorption cross-sections of ozone in the ultraviolet and




Figure 9. Residual of the quadratic fit




Figure 10. Residual of the cubic fit

Table 1. Emission lines of the Discharge Lamps used for Brewer calibration

Lamp	Line (nm)	Slits
Mercury (Hg)	289.36	01
Hg	296.728	03
Zinc (Zn)	301.836	05
Zn	303.578	05
Cd (multiplet)	308.082	0-5
Cd	313.3167	05
Cd	326.1055	05
Zn	328.233	05
Hg	334.148	05
Cd	340.3652	05
Cd	349.995	45
Cd (multiplet)	361.163	5

Table 2. Ozone Absorption Coefficient in atm cm^{-1} calculated using four absorption cross sections,

	Parametrized	Measured
BRW	0.3381	0.3407
B&P	0.3330	0.3360
DMB	0.3483	0.3514
SDK	0.3392	0.3422

Weber, M., Gorshelev, V., and Serdyuchenko, A.: Uncertainty budgets of major ozone absorption cross sections used in UV remote sensing applications, *Atmospheric Measurement Techniques*, 9, 4459–4470, <https://doi.org/10.5194/amt-9-4459-2016>, <https://www.atmos-meas-tech.net/9/4459/2016/>, 2016.

visible spectral regions: Status report 2015, *Journal of Molecular Spectroscopy*, 327, 105–121, 2016.

Savastiouk, V.: Improvements to the direct-sun ozone observations taken with the Brewer spectrophotometer, 2005.

⁵ Serdyuchenko, A., Gorshelev, V., Weber, M., Chehade, W., and Burrows, J. P.: High spectral resolution ozone absorption cross-sections: Part II. Temperature dependence, *J. Quant. Spectrosc. Radiat. Transfer*, 2012.

Table 3. Ozone Absorption Coefficient in atm cm^{-1} calculated using four absorption cross sections,

	brw_scan	ptb_scan	opo_quad	opo_cubic	lamp_quad	lamp_cubic
SGW	0.3409	0.3368	0.3442	0.342	0.3446	0.3412
ratio	1	0.9881	1.0096	1.0033	1.0108	1.001 height